

ANY TIME ANY PLACE

All-in-one kit for intermittent catheterization.
Easy for beginners. Great for long-term use.

LoFric® Hydro-Kit™

A man with a beard and short dark hair is sitting in a wheelchair on the deck of a boat. He is wearing a dark blue patterned button-down shirt with the sleeves rolled up and grey cargo pants with the cuffs rolled up. He is smiling and looking out towards the right side of the frame. The background shows the rigging of several sailboats under a clear blue sky. The text is overlaid on the right side of the image.

“Self-catheterization
shouldn’t stop you
from doing what
you love to do.”

CATHETERIZATION WHEN AND WHERE YOU WANT

Having to use a catheter certainly shouldn't stop you from doing what you want to do. And with a catheter kit, you have everything you need for emptying your bladder without access to a bathroom. A catheter and a collection bag, all in one kit.

Catheter kits are useful both when you are out and about, and when you have difficulty leaving your bed or wheelchair.

Discreet. No fuss. No touch.

At Wellspect, we constantly strive to improve our products to meet our users' needs. We've worked with catheters, in close consultation with users, for over 30 years. That way, we have learned what really matters.

Our users have taught us that a catheter kit should be easy to use with a minimum of fuss. It should have a non-touch function to avoid bacteria and possible infections. And it should be discreet - small enough to carry and easy to dispose of.



LoFric Hydro-Kit is the answer

LoFric® Hydro-Kit™ is a complete all-in-one catheter kit designed to address both your medical and lifestyle needs. It is easy to bring and use anywhere - without ever making concessions to health or safety.



A COMPLETE ALL-IN-ONE CATHETER KIT

The single use catheter kit that combines a hydrophilic catheter, salt solution sachet and collection bag.



NEW Refined grip


Textured Insertion Grip ensures a firm hold and promotes hygiene thanks to a non-touch technique.

The loops make it easy to use, even with poor hand function.



LoFric catheter with Urotonic™ Surface Technology

This unique coating makes the hydrophilic catheter surface smooth, comfortable and safe during both insertion and withdrawal.







Collection bag
With printed scale (1000 ml).

Opening for emptying the collection bag

Loop for hanging or holding the collection bag



NEW

Instant activation

Burst the sachet with salt solution, soak the catheter and it is immediately ready to use.

Long neck

Allows you to hang the catheter over your knee for convenient catheterization.

Lets you seal the collection bag easily after use, simply by tying a knot.

WHAT USERS LIKE ABOUT LOFRIC HYDRO-KIT



“Small enough to put in my pocket”

“I live an active life, travelling a lot, and sometimes I find myself in places without accessible toilets. That is when I appreciate using LoFric Hydro-Kit. It’s small enough to put in my pocket, and thanks to the built-in collection bag I can more or less use any private room to empty my bladder. Afterwards I put it in my backpack for disposal whenever suitable. I prefer a kit over separate components that must be assembled, especially when I’m travelling as it greatly reduces the risk of leakage.”

Kent, paraplegic



“It simplifies my life”

“What I love about LoFric Hydro-Kit is that it simplifies my life. I can empty my bladder directly from my wheelchair and sometimes even in bed. And it is a kit, all in one. My hands are not very strong, and the loops on the product help me both to open and handle the product.”

Mona, tetraplegic

GET IT RIGHT FROM THE START

Complications like urethral trauma, strictures and urinary tract infections (UTIs) might not be noticeable from the outset. Therefore, to ensure long-term urethral health, make sure to choose the right catheter – right from the start.

Safe for long-term use

LoFric catheters are coated with our unique Urotonic™ Surface Technology, which means that LoFric's surface has the same properties as your urethra – the salt concentration is the same as in urine. This ensures that the catheter remains smooth, comfortable and safe, during both insertion and withdrawal.



Scientifically proven

LoFric catheters are proven in multiple studies to be safe, also for long-term use.



0% trauma after 5-9 years with LoFric.¹



10-20 times lower expected risk of strictures.²



64% lower expected risk of UTIs.³



48% lower expected risk of hematuria.⁴



Learn more about LoFric products

References

1. 7 years follow up of 30 LoFric patients. Waller et al. J Urol 1995;153:345-8
2. Survey in 288 long-term (10 years) LoFric users. Hakansson et al. Urol Nurs. 2015;35:239-47
3. Meta-analysis of 464 patients. Li et al. Arch Phys Med Rehabil. 2013;94:782-7
4. Health Technology Assessment (HTA). Rognoni and Tarricone, CERGAS Bocconi University 2016

At Wellspect we value the people behind our success as a leading provider of life-changing products for bowel and bladder management such as the well-known brands LoFric® and Navina™. From the thousands of users and healthcare professionals worldwide who inspire our innovative solutions, we know that working together is the best way to advance continence care, giving our users more time for life. Building on over 30 years of life-improving performance, we passionately strive to make a difference every day to everyone who needs our products and services.

Wellspect. A Real Difference.

wellspect.com

