

Navina Mini Product Information

Product description

Low volume (185ml) transanal Irrigation device.

Shelf life and uses

All Navina Mini components have a shelf life of 2 years. The cone is a single use product and should only be used once. The water container and optional extension tube can be used 60 times.

Material composition

Main product components: water container, rectal cone and extension tube are made of thermoplastic elastomer (TPE).

Cone Hydrophilic Surface Coating: Polyvinylpyrrolidone (PVP).

Primary Packaging Material: Polyethylene (PE)

Secondary Packaging Material: Cardboard and Corrugated cardboard.

Navina Mini **are not manufactured** with latex, PVC or phthalates.

Intended use

Navina Mini is intended for Transanal Irrigation by instilling a low volume of water into rectum through a single use rectal cone.

Indications for use

Navina Mini is indicated for adults, adolescents and children from 3 years who suffer from incomplete emptying and/or fecal incontinence.

Contraindications

Do NOT use Navina Mini if you have one or more of the following:

- Known anal stenosis.
- Active inflammatory bowel disease.
- Acute diverticulitis and diverticular abscess.
- Colorectal cancer.
- Ischemic colitis.
- You are within three months of anal or colorectal surgery.

Since the list may not be exhaustive, healthcare professionals should always consider individual user factors as well.



Navina Mini

Reference	Type	Specification
6931000	Navina Mini Start Set	1 water container, 5 rectal cones, 1 extension tube
6930700	Navina Mini Set	1 water container, 30 rectal cones
6930800	Navina Mini Extension Set	1 water container, 30 rectal cones, 1 extension tube
6930900	Navina Mini Refill Set	30 rectal cones

Navina™ Mini



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00. www.wellspect.com

CE 2797

Rx Only