

# LoFric Elle Pro catheters

## – Product data sheet

**LoFric® Elle™ Pro**, a single use hydrophilic, ready to use catheter with multiple eyelets, designed to allow for complete bladder emptying without the need for repositioning.

### Sterilization and shelf life

Sterilized by irradiation. Use before expiry date on the package.

### Package

30 catheters per customer box.  
2 customer boxes per shipping box.  
5 catheters per sample pack.

### Material composition

Hydrophilic surface coating: Polyvinyl pyrrolidone (PVP).

Catheter: Polyolefin-based elastomer (POBE).

Wetting fluid: Water-based wetting fluid with glycerol.

Item packaging: Polypropylene.

### Product range

Reference	Type	Length	Size	No. eyelets
419 10 XX*	LoFric® Elle™ Pro Nelaton (female)	10 cm**	CH 10	12
419 12 XX*	LoFric® Elle™ Pro Nelaton (female)	10 cm**	CH 12	12
419 14 XX*	LoFric® Elle™ Pro Nelaton (female)	10 cm**	CH 14	12

\*Generic article number without the two digit suffix specific for a region or country destination when distributing an article.

\*\*Effective length, length of the catheter that can be inserted into the body.

All products available in sample pack, 5 catheters. To order, replace the last two digits in the Reference number with 05.

## LoFric Elle Pro



- **LoFric Elle Pro multiple eyelets** are designed to allow for complete bladder emptying without the need to reposition the catheter.<sup>1</sup>
- **Improving the success rates** of catheterization – the L-shaped handle improves visibility when using a mirror and makes it easier to locate the urethra. It also supports a more upright, natural position during catheterization.

- **L-shaped handle** for a good grip, extra reach, and non-touch technique.
- **Ready to use** – ready when you need it, to enable fast handling and catheterization.
- **Smooth eyelets.**
- **LoFric catheters** are coated with the unique Urotonic™ Surface Technology. This ensures

the catheter remains smooth and safe, both during insertion and withdrawal.<sup>2,3</sup>

- **Lower container** is recyclable as plastic (polypropylene).
- **Not manufactured with PVC**, phthalates or latex.



**LoFric® Elle™ Pro**

## References

1. Data on file.
2. Waller et al., Spinal Cord, 1997.
3. Dunger et al., Acta Endocrinologica, 1988.

At Wellspect, we develop innovative products and solutions that make a real difference in the lives of people managing continence care needs, empowering them to live each day to the fullest. We have been a trusted leader in healthcare for over 40 years with our product brands LoFric® and Navina™, and we continue to expand access to our user-friendly solutions. Every day, we strive to maintain the highest quality standards while minimizing the environmental impact of our products. By working closely with users and healthcare professionals, we aim to improve clinical outcomes and enhance quality of life for people around the world.

**Wellspect. A Real Difference.**

For more information about our products and services, please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on LinkedIn, X, Facebook and Instagram.

[wellspect.com](https://www.wellspect.com)

Wellspect AB, Aminogatan 1, SE-431 21 Mölndal, Sweden. Phone: +46 (0)31 376 40 00.

LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.  
Phone: +46 31 376 40 00. [www.wellspect.com](https://www.wellspect.com)

CE 2797

R<sub>x</sub> Only