

**LoFric® Origo™ Pro**, Single use, hydrophilic urinary catheter, with salt solution for activation and multiple eyelets, designed to allow for complete bladder emptying without the need for repositioning.

#### Sterilization and shelf life

Sterilized by irradiation. Use before expiry date on the package.

#### Package

30 catheters per customer box.  
4 customer boxes per shipping box.  
5 catheters per sample pack.

#### Material composition

Hydrophilic surface coating:  
Polyvinyl pyrrolidone (PVP).

Catheter: Polyolefin-based elastomer (POBE).

Sachet: Laminate foil of polyethylene terephthalate (PET), aluminum and polyethylene (PE). Sterile water with sodium chloride (NaCl).

Sleeve: Polyethylene (PE)

#### Item packaging:

Top foil: Laminate foil of polyethylene terephthalate (PET) and polyethylene (PE).

Bottom foil: Polyolefin.

#### Product range


| Reference  | Type                 | Length | Size  | No. eyelets |
|------------|----------------------|--------|-------|-------------|
| 441 12 XX* | LoFric® Origo™ Pro** | 40 cm  | CH 12 | 12          |
| 441 14 XX* | LoFric® Origo™ Pro** | 40 cm  | CH 14 | 12          |
| 441 16 XX* | LoFric® Origo™ Pro** | 40 cm  | CH 16 | 12          |

\*Generic article number without the two digit suffix specific for a region or country destination when distributing an article.

\*\* LoFric Origo Pro contains Sleeve, Flexible tip and multiple eyelets.

All products available in sample pack, 5 catheters. To order, replace the last two digits in the Reference number with 05.

## LoFric Origo Pro



- **LoFric Origo Pro multiple eyelets** are designed to allow for complete bladder emptying without the need to reposition the catheter.<sup>1</sup>
- **Loop opening label.**
- **Adhesive tab** on back for reliable hanging.
- **Instant activation** – sterile salt solution included, squeeze and it's ready to use.
- **Protective sleeve** for excellent control allows for non-touch technique.
- **Smooth eyelets.**
- **LoFric Origo catheters** are coated with the Urotonic Surface Technology. The catheter remains smooth, both during insertion and withdrawal.<sup>2,3,4</sup>
- **Not manufactured with PVC,** phthalates or latex.



## References

1. Data on file at info@wellspect.com.
2. Waller et al., 1997.
3. Lundgren et al., 2000;
4. Johansson et al., 2013

At Wellspect, we develop innovative products and solutions that make a real difference in the lives of people managing continence care needs, empowering them to live each day to the fullest. We have been a trusted leader in healthcare for over 40 years with our product brands LoFric®, Navina™ and Surity™, continuously expanding access to our user-friendly solutions. Every day, we strive to maintain the highest quality standards while minimizing the environmental impact of our products. By working closely with healthcare professionals and people who rely on our solutions, we aim to improve clinical outcomes and enhance quality of life for people around the world.

**Wellspect. A Real Difference.**

For more information about our products and services, please visit [Wellspect.com](https://wellspect.com).

Join the conversation on LinkedIn, X, Facebook and Instagram.

**[wellspect.com](https://wellspect.com)**

Wellspect AB, Aminogatan 1, SE-431 21 Mölndal, Sweden. Phone: +46 (0)31 376 40 00.

**LoFric**



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.  
Phone: +46 31 376 40 00. [www.wellspect.com](https://www.wellspect.com)

**CE** 2797

**R<sub>x</sub>** Only